

## L3-SERIES PIPETTES



### PIPETTES-L

Continuously adjustable pipette with tip ejector designed to measure and transfer liquids precisely and safely.

10 models in optimal volume ranges permit accurate measuring of any volume between 0,1 $\mu$ l to 10ml.

#### Characteristic Features:

- Perfect results of measurements
- Excellent ergonomical shape of the handle makes the work easy and pleasant
- Easy to read digital display
- Easy set up tip ejector's length for any kind of desired tip
- Easy removal of tip ejector when working with narrow tubes
- May be re-calibrated for solutions of any density and viscosity
- Each unit comes with an individual quality certificate, holder and re-calibration key
- Autoclavable shaft
- Maintenance free
- The pipette and its components are made of the best plastic and the best stainless steel and can be used safely with acid and alkaline solution

| Cat. No. | Model   | Volume Range $\mu$ l | Accuracy error % | Precision error % |
|----------|---------|----------------------|------------------|-------------------|
| 3.1.1.   | L2      | 0,1-2                | 12 - 1,5         | 6,0 - 0,7         |
| 3.1.2.   | L10     | 0,5-10               | 4 - 0,5          | 4,0 - 0,4         |
| 3.1.3.   | L20     | 2-20                 | 4 - 0,8          | 3,0 - 0,4         |
| 3.1.3.0. | L50     | 5-50                 | 3 - 0,8          | 3,0 - 0,4         |
| 3.1.4.   | L100    | 10-100               | 1,6 - 0,8        | 0,3 - 0,3         |
| 3.1.5.   | L200    | 20-200               | 1,2 - 0,6        | 0,6 - 0,3         |
| 3.1.5.0. | L250    | 50-250               | 1,2 - 0,6        | 0,6 - 0,3         |
| 3.1.6.   | L1.000  | 100-1000             | 0,9 - 0,6        | 0,5 - 0,25        |
| 3.1.7.   | L5.000  | 1000-5000            | 0,6 - 0,5        | 0,3 - 0,2         |
| 3.1.8.   | L10.000 | 1000-10000           | 0,5 - 0,4        | 0,2 - 0,1         |